



# PRODUCT SPECIFICATION



**Product:** CEDERROTH FIRST AID & BURN STATION

**Reference code:** 490960



## 1. General description

A First Aid station containing a variety of first aid products. All products are easily available and protected against dust and moisture by transparent lids. Space for additional products is provided behind a lid with first aid instructions. Adhesive bandages (plasters) and wound cleansers are placed in dispensers for easy access and theft-proof storage. A special key is provided for refill replacement.

## 2. Intended use

To be used in first aid situations, especially suitable for work places where there is a risk of burns.

## 3. Technical description

### *Product dimensions*

(l x w x h) cm: 56 x 29 x 12

### *Material*

First aid station	Green parts:	Polystyrene, ABS-plastic
	Transparent lid:	Polycarbonate
	Sealing profiles:	EPDM rubber (synthetic rubber)
First aid products		Refer to Product specifications for each product

## 4. Composition

### *Consumer package*

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>
51616965	First Aid station (empty)	1
1910	Cederroth 4-in-1	2
1911	Cederroth 4-in-1 Mini	3
901900	Burn Gel Dressing 2-p	2
901901	Burn Gel Spray	2
6735CAP	Salvequick Blue Detectable Plasters, regular, in dispenser	1
6754CAP	Salvequick Blue Detectable Plasters, regular/fingertip, in dispenser	1
676100	Cederroth Soft Foam Bandage, blue, in dispenser	1
323700	Salvequick Wound Cleanser 20-p	1
1921	Cederroth Breathing mask	1
51616417	First Aid Instruction (label)	1
51616620	Key to the dispensers	1
51616668	Refill instruction with drill template	1



# PRODUCT SPECIFICATION



## *Transport package*

Number of consumer packages in transport package:	1
Dimensions (l x w x h) cm:	57 x 29 x 12
Weight (kg):	3,3

## 5. Labeling

### *On product*

Lot number and serial number  
Instructions  
Reference to the products

### *On consumer/transport package*

Cederroth address, telephone- and  
faxnumber, web-address  
Product name, Reference code  
Lot-number, Expiry date, EAN-code  
Quantity.

## 6. Storage conditions

The products in the first aid station should be stored in their original packages for maximum protection until time of use.

Shelf life: 5 years from production date for each component

## 7. Accessories

No accessories are necessary for the use of this product.

## Authorization: