

Product: CEDERROTH EYE WASH
5x500 mL IN A PORTABLE CASE

Reference code: 7255



1. General description

A plastic case with transparent lid and a handle.
Contains 5 bottles of Cederroth Eye Wash, 500 mL.

2. Intended use

The Cederroth Eye Wash is intended for washing away foreign particles and dangerous liquids, like alkali or acids, from an eye in an emergency first aid situation at work places. The Cederroth buffered Eye Wash also has a neutralising effect on splashes from alkalis and acids. The effect is in general stronger for alkalis than for acids.

The Cederroth Eye Wash is intended for the first acute phase at the very accident site before getting to a plumbed or self-contained eye wash equipment or seeking medical attention.

The portable case is intended for mobile work places and for continuous irrigation during transport to the hospital.

3. Technical description

Product dimensions : l x w x h (cm) 40 x 28 x 8,5 plus handle

Material

Plastic case	Lid : Transparent ABS plastic Box: Green Polystyrene Insert: White Polystyrene
Eye wash	Refer to product specification for REF 7251

4. Sterilization

Sterilization method: Aseptic filling of bottles

5. Composition

Consumer package

Transport package

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>	Number of consumer packages
7251	Cederroth Eye Wash	5	in transport package: 3
51618723	Plastic case	1	Dimensions (l x w x h) cm: 41 x 27 x 29
29164200	Insert	1	Weight (kg): 12,6

6. Labeling

On product/ consumer package

Cederroth address, telephone- and fax-number, e-mail and web-address
Product name, Reference code
Lot-number, Expiry date, EAN-code
Instructions for use
CE +NB#, Storage conditions

On transport package

Cederroth address, telephone- and faxnumber, web-address
Product name, Reference code
Lot-number, Expiry date, EAN-code
Quantity
Symbols: LOT, USE BY

7. Storage conditions

Store in temperatures between +2°C and +40°C. Do not freeze. Store in the original packaging for maximum protection until time of use.

Shelf life: 4,5 years from production date (eye wash)

8. Accessories

The following accessories can be used with the product but are not necessary.

190400

Wall bracket for First Aid Cases



Authorization: