



# PRODUCT SPECIFICATION



**Product:** SAVETT WOUND CLEANSER, REFILL

**Reference code:** 3227



## 1. General description

Sterilized sachets containing non-woven wipes saturated with wound cleansing liquid. The wipes are individually packed in sealed, laminated foil, which prevents leakage and guarantees sterility. The sachets are packed in a refill box, intended for a Savett Wound cleanser dispenser.

## 2. Intended use

The Savett Wound Cleanser is intended for mechanical cleansing of minor wounds when clean tap water is not available.

## 3. Technical description

| <i>Product dimensions</i> | <i>(l x w x h) mm</i> | <i>Material</i> |                                  |
|---------------------------|-----------------------|-----------------|----------------------------------|
| Wipe :                    | 120 x 110             | Wipe:           | White non-woven tissue (Norafin) |
| Sachet                    | 60 x 58               | Sachet (foil):  | Paper/Aluminium/Surlyn           |
| Box:                      | 60 x 60 x 151         | Liquid:         | Physiological saline solution    |

## 4. Sterilization

Sterilization method: Irradiation

## 5. Composition

### *Consumer package*

| <u>Ref#</u> | <u>Product/Material</u> | <u>Pcs</u> |
|-------------|-------------------------|------------|
| 3239        | Wound cleanser sachet   | 40         |

### *Transport package*

|   |              |
|---|--------------|
| Number of consumer packages in transport package: | 10           |
| Dimensions (l x w x h) cm:                        | 33 x 14 x 16 |
| Weight (kg):                                      | 1,9          |

## 6. Labeling

### *On product/ consumer package*

Cederroth contact information, Product name, Reference code. Lot-number, Expiry date, EAN-code. Do not reuse, CE 0413, Sterile R

### *On transport package*

Cederroth contact information, Product name, Reference code, Lot-number, Expiry date, EAN-code, Quantity.

## 7. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

Shelf life: 5 years from production date

## Authorization: